

PRODUCT SPECIFICATIONS

Product Name:	Human CA 15-3 Standard Stock Solution
Catalog Number:	90008
Source:	Human BTA cell line supernatant
Lot Number:	
CA 15-3 Concentration:	
Quantity:	
Expiration Date:	
Product Form:	In buffer solution containing 1.0% (v/v) ProClin-300 as preservative
Storage:	-20°C or below.
Viral Tests:	The source materials used for manufacture of this product tested negative for HBsAg, HIV 1&2 and HCV by FDA-approved methods. No methods can completely assure the absence of these agents and other infectious agents may be present.