

Product Name: Mouse Anti-human CA153 monoclonal antibody , clone 12A9
Product: Purified Monoclonal Antibody **Catalog No.:** DM0005-B5
Lot Number: **Concentration:**
Production Date:
Exp Date:

Antibody Purification

Starting Material: cell culture supernatant
Axial Flow Column: Protein G Sepharose
Buffer/Components: PBS (pH7.3), 0.1%NaN3
Storage Conditions: Store at 2-8°C ,Do Not freeze.

Test:	Method:	Document:	Results:
Protein Concentration	BCA	G3576	
Purity	SDS-Page (reduced)	G3302	
Isotype	ISOSTrip™	G3346	
Application	ELISA	Datasheet	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date