

Product Name: SCC protein**Catalog No.:** DP0035-P1**Lot Number:****Concentration:****Manufactured:****EXP date:****Antibody Purification****Buffer/Components:** PBS(PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide**Storage Conditions:** -20°C. Avoid repeated freeze-thaws.

Test:	Method:	Document:	Results:
Protein Concentration	A280	G3476	
Purity	SDS-Page (reduced)	G3302	
Application	ELISA	Datasheet	

*The data from the raw material without adding the dilution buffer

Prepared By:_____**Quality Control-Date****Reviewed By:**_____**Quality Assurance-Date**

No.168 Mashan Meiliang Road | Wuxi, Jiangsu Province, China |0510-81830666

www.origene.com.cn techsupport@origene.cn