

Product Name: CHD4(Mi-2beta)
Product: Purified Recombinant Protein
Lot Number:
Production Date:

Catalog No.: DP0097-P1
Concentration:
Expiring Date:

Protein Purification

Starting Material: Cell Lysate
Buffer/Components: 25 mM Tris-HCl, 100 mM glycine, 250 mM NaCl, 10% glycerol, pH 7.3
Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date