

Product Name: AChR (CHRNA3)
Product: Purified Recombinant Protein
Lot Number:
Production Date:

Catalog No.: DP0114-P1
Concentration:
Expiring Date:

Protein Purification

Starting Material: cell culture supernatant
Buffer/Components: PBS, 10% glycerol
Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date