

Product Name: MORC3
Product: Purified Recombinant Protein
Lot Number:
Production Date:

Catalog No.: DP0118-P1
Concentration:
Expiring Date:

Protein Purification

Starting Material: Cell Lysate
Buffer/Components: 25 mM Tris-HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage Conditions: Store at -80°C. Avoid repeated freeze-thaw cycles.

Test:	Method:	Document:	Results:
Protein Concentration	Bradford protein assay method		
Purity	SDS-Page (reduced)	G3302	

Prepared By: _____
Quality Control-Date

Reviewed By: _____
Quality Assurance-Date

样品信息																								浓度														纯度														WB																																
序号	编号	送检组别	送检日期	样品TP号 (SKU)	批次号	对应项目号	样品类型	Symbol	流水号	RC号	理论分子量kDa	生产体系	蛋白定位	研发or生产	项目归属	优先程度	送检体积ul	蛋白浓度检测方法	空白buffer	SDS-PAGE建议上样量	Western Blot建议上样量	蛋白标签	备注	送检人员	QC接收人	QC检测日期	检测次数	是否24h内检测	24h内未安排原因	稀释倍数	QC实验人	QC检测方法	检测浓度(ug/ul)	QC反馈人	QC反馈日期	QC备注	QA反馈人	QA反馈日期	QC检测日期	检测次数	是否24h内检测	24h内未安排原因	QC检测人	稀释方法	处理方法	上样体积ul	BSA对照批次	工作单号	工作单链接	泳道	QC反馈人	QC反馈日期	QC备注	初判结果	初判纯度分析	初判结论	初判人	初判反馈日期	初判反馈人	QA复核人	复核结果	复核纯度分析	复核结论	QA反馈日期	QC检测人	稀释方法	上样体积ul	孵育抗体TA-SKU	孵育抗体批次	工作单号	工作单链接	泳道	QC反馈人	QC反馈日期	QC备注	初判结果	初判结论	初判人	初判反馈日期	初判反馈人	QA复核人	复核结果	复核结论	QA反馈日期
P1-W	1	白部真核生产	2025/4/2	DP0118-P1	WX03D222	8-P1/TP7	蛋白	MORC3	CA93	CA93	56.1	真核	胞内蛋白	自免备货	研发		70ul	G250	空白1: 25m	5	2	FLAG		焦灿	刘莉	2025/4/2	1	是		1	王紫洁	G250	0.21	刘莉	2025/4/2	检样品均为7ul/thermo		2025/4/2	1	是		刘莉/王紫洁	原液+7.5ul5min+1W转离	5	20250402	20250402-1	20250402-1	4	王紫洁	2025/4/2		主带位置59k	80.4%	合格	苏晓烽	2025/4/3	苏晓烽	李欣佩	主带位置60k	81.1%	合格	2025/4/3	仲浩成	原液+5ul15min	1	TA180144	F002	B20250403	B20250403-1	2	仲浩成	2025/4/3		主带位置59k	合格	苏晓烽	2025/4/7	苏晓烽	李欣佩	主带位置60k	合格	2025/4/7		